

MAT KO Transformation

Introduction

ste4 Δ ::TRP1 yeast strain derived from SEY α or SEY α are transformed with [MAT KO PCR product](#) to replace MAT α or MAT α with His3. A heat shock time of 1 h increases transformation efficiency for SEY derived strains when compared to 20 min heat shock.

Materials

- > YPAD
- > 50% PEG
- > 1 M LiAc
- > 2mg/ml ssDNA
- > [MAT KO PCR product](#) (at least 30 ng/ μ L)
- > Sterile H₂O
- > Water bath or heating block at 42°C
- > Disposable 17x100 mm culture tubes with closures (VWR Catalog # 60818-725)
- > [SC-His-Trp plates](#)

Procedure

Competent Cells

Day 1

1. Overnight culture (ONC): Inoculate 2 mL of YPAD with a single yeast colony and incubate culture on the shaker at 30°C/200 rpm/overnight.

Day 2

2. Inoculate 2.5 mL of YPAD in disposable 10 mL culture tubes with 200 μ L of ONC and grow it on the shaker at 30°C/200 rpm/4 h.

04:00:00



3. Set water bath or heating block to 42°C.
4. Indicate number of transformations and PCR product concentration:

Settings		
	# of Transformations	MAT KO PCR product [ng/ μ L]
1	24	100



5. If the carrier DNA has not been denatured yet, denature the carrier DNA in heating block at 95°C for 5 min and chill immediately on ice before use. If carrier DNA has already been denatured then thaw it on ice.

6. Prepare 0.1 M LiAc, freshly diluted from 1 M LiAc:

0.1 M Li Ac		
	Materials	mL
1		26 X [mL]
2	1 M LiAc	2.6
3	ddH ₂ O	23.4
4	Total	26



7. Centrifuge 3,000 g/5 min; discard supernatant and resuspend cell pellet in 1 mL sterile H₂O.

8. Centrifuge 3,000 g/5 min; discard supernatant and resuspend cells in 1 mL of 0.1 M LiAc.

9. Centrifuge 3,000 g/5 min; discard supernatant.

Transformation

10. Prepare transformation master mix by adding reagents in order (pipetting 50% PEG slowly), two additional reactions are included to account for pipetting discrepancies:

Transformation master mix		
	A	B
1		26 X [μ L]
2	50% PEG	6240
3	1 M LiAc	936
4	ssDNA (2 mg/ml)	1300
5	MAT KO PCR product	260
6	ddH ₂ O	624
7	Total	9360



11. Add 360 μ L of transformation master mix to each tube and immediately vortex until the cells are smoothly suspended.

12. Heat shock cells in water-bath or heat-block: 42°C/1 h.



13. Centrifuge 3,000 g/5 min; discard supernatant and resuspend each cell pellet in 100 µL sterile H₂O, then plate cells on the SC-His-Trp.

Additional resources

PCR product amou...		^
	KO PCR prod- uct [ng]	
1	1000	

Transformation single reaction			^
	A	B	
1		[µL]	
2	50% PEG	240	
3	1 M LiAc	36	
4	2mg/ml ssDNA	50	
5	PCR product	10.0	
6	dHd2O	24.0	
7	Total	360	