

ste4 Colony PCR

Introduction

For genotyping of ste4 locus: the primers bind upstream (F) and downstream (R) of ste4, and PCR is used to determine if TRP1 has replaced ste4 by homologous recombination.

Materials

- > TaKaRa LA Taq® DNA Polymerase (Mg2+ free buffer) RR002A
- > 10 µM Primer F: OLPr022 STE4 -133 F
- > 10 µM Primer R: OLPr023 STE4 1402 R
- > [SC-Trp plates](#)
- > Template: yeast colonies transformed with [ste4 KO PCR product](#)

Procedure

PCR

1. Chose the number of 10 µL reactions (include 3 additional reactions to account for wild-type, mutant, and no-template control):

Settings	
	# of reactions
1	0.5

2. Combine the following reagents to prepare the master mix (contains 0.5 additional reaction to account for pipetting discrepancies), and aliquot 10 µL master mix/pcr tube:

Master mix		
	A	B
1		1 X [µL]
2	ddH2O	5.9
3	10X LA PCR Buffer	1
4	25 mM MgCl2	1
5	dNTP Mix (2.5 mM each)	1.6
6	10 µM OLPr022 STE4 -133 F	0.2
7	10 µM OLPr023 STE4 1402 R	0.2
8	LA Taq Polymerase (5 U/µL)	0.1
9	Total	10

3. In sterile work environment (e.g. laminar flow), gently touch a yeast colony with a 10 µL pipette tip so a small number of cells stick to it, touch the tip to a fresh reference SC-Trp plate, and resuspend cells in PCR reaction.

Do not use too many cells as they might inhibit the PCR reaction. The PCR reaction should only look slightly opaque from the resuspended cells. If the reaction looks milky then you might have used too many cells.

The reference plate is incubated 30°C/2–3 days.

4. PCR program "STE4COLO" in folder "STEPHAN":

94°C/5 min—[94°C/0.5 min—61°C/0.5 min—68°C/1.5 min]x35—68°C/5 min—10°C/∞

Gel

5. Gel: 10 µL/0.5 µg O'GeneRuler DNA Ladder Mix/1% TAE/120 V/35 min

00:35:00



Expected product sizes:

ste4 1535 bp

ste4Δ::TRP1 1232 bp

Stop processing any yeast strains with the wrong genotype (i.e. discard any corresponding overnight cultures and patches on reference plate).

Additional resources

Primers and template		
	A	B
1	Forward Primer	OLPr022 STE4 -133 F
2	Reverse Primer	OLPr023 STE4 1402 R
3	Template	yeast colonies



PCR Volume and reagent concentrations

^

	A	B
1	PCR Volume [μ L]	10
2	Starting buffer conc [X]	10
3	Final buffer conc [X]	1
4	Starting Mg ²⁺ conc [mM]	25
5	Final Mg ²⁺ conc [mM]	2.5
6	Starting dNTP conc each [mM]	2.5
7	Final dNTP conc each [mM]	0.4
8	Starting Primer F conc [μ M]	10
9	Final Primer F conc [μ M]	0.2
10	Starting Primer R conc [μ M]	10
11	Final Primer R conc [μ M]	0.2
12	Starting Polymerase conc [U/ μ L]	5
13	Final Polymerase conc [U/ μ L]	0.05
14	Starting Template conc [ng/ μ L]	5
15	Final Template conc [ng/ μ L]	0

Single PCR reaction

^

	A	B
1		[μ L]
2	ddH ₂ O	5.9
3	10X PCR Buffer	1
4	25 mM MgCl ₂	1
5	dNTP Mix (2.5 mM each)	1.6
6	10 μ M Primer F	0.2
7	10 μ M Primer R	0.2
8	LA Taq Polymerase (5 U/ μ L)	0.1
9	Template (5 ng/ μ L)	0
10	Total	10