

Plasmid Transformation

Introduction

MAT Δ ::His3 STE4 Δ ::TRP1 yeast strain derived from SEY α or SEY α are transformed with OLP003 or OLP004. A heat shock time of 1 h increases transformation efficiency for SEY derived strains when compared to 20 min heat shock.

Materials

- > YPAD
- > 50% PEG
- > 1 M LiAc
- > 2mg/ml ssDNA
- > OLP003
- > OLP004
- > Sterile H₂O
- > Water bath or heating block at 42°C
- > Disposable 17x100 mm culture tubes with closures (VWR Catalog 60818-725)
- > SC-Leu plates
- > SC-Ura plates

Procedure

Competent cells

Day 1

1. Overnight culture (ONC): Inoculate 2 mL of YPAD with a single yeast colony and incubate culture on the shaker at 30°C/200 rpm/overnight.

Day 2

2. Inoculate 5 mL of YPAD in disposable 10 mL culture tubes with 400 μ L of ONC and grow it on the shaker at 30°C/200 rpm/4 h

This will make 2 x 1 mL competent cells.

3. Set water bath or heating block to 42°C.
4. Indicate number of transformations and plasmid concentration:

Settings				^
	# of Transformations	OLP003 [ng/ μ L]	OLP004 [ng/ μ L]	
1	24	610	550	

5. If the carrier DNA has not been denatured yet, denature the carrier DNA in heating block at 95°C for 5 min and chill immediately on ice before use. If carrier DNA has already been denatured then thaw it on ice.

6. Prepare 0.1 M LiAc, freshly diluted from 1 M LiAc:

0.1 M Li Ac			^
	Materials	mL	
1		26 X [mL]	
2	1 M LiAc	5.2	
3	sterile H ₂ O	46.8	
4	Total	52	

7. Centrifuge 3,000 g/5 min; discard supernatant and resuspend cell pellet in 2 mL sterile H₂O.

8. Centrifuge 3,000 g/5 min. Discard supernatant and resuspend cells in 2 mL of 0.1 M LiAc. Split into 2 x 1 mL in microcentrifuge tubes (transform one with OLP003 and one with OLP004 below).

9. Centrifuge 3,000 g/5 min (use 7,500 rpm in microcentrifuge); discard supernatant.

10. Prepare one transformation master mix for each plasmid by adding reagents in order (pipetting 50% PEG slowly), two additional reactions are included to account for pipetting discrepancies:

Transformation master mix				^
	A	C	C	
1		26 X OLP003 [μ L]	26 X OLP004 [μ L]	
2	50% PEG	6240	6240	
3	1 M LiAc	936	936	
4	ssDNA (2 mg/ml)	1300	1300	
5	Plasmid	42.6	47.3	
6	ddH ₂ O	841.4	836.7	
7	Total	9360	9360	

11. Add 360 μ L of transformation master mix to each tube.

12. Vortex tubes until the cells are smoothly suspended.

13. Heat shock cells in water-bath or heat-block: 42°C/1 h.



14. Centrifuge 3,000 g/5 min (use 7,500 rpm in microcentrifuge); discard supernatant and resuspend each cell pellet in 100 μ L sterile H₂O, then plate cells on SC-Leu for OLP003 transformations or SC-Ura for OLP004 transformations.

Additional resources

Plasmid amount pe... [^]	
	Plasmid [ng]
1	1000

Transformation single reaction [^]			
	A	B	C
1		OLP003 [μ L]	OLP004 [μ L]
2	50% PEG	240	240
3	1 M LiAc	36	36
4	2mg/ml ssDNA	50	50
5	PCR product	1.6	1.8
6	dHd2O	32.4	32.2
7	Total	360	360