

MAT KO PCR

Introduction

The primers contain homologous sequences upstream (F) and downstream (R) of MAT and bind to the Universal F1 or Universal R1 region flanking His3 on OLP45. The PCR product is used to delete MAT locus with His3.

Materials

- › TaKaRa LA Taq® DNA Polymerase (Mg2+ free buffer) RR002A
- › 10 µM OLPr008 MAT -33 UKO F
- › 10 µM OLPr012 MAT +34 SEY UKO R
- › Template: [OLP045](#)

Procedure

PCR

- ✓ 1. Chose the number of 50 µL reactions (include one additional reaction for no-template control):

Settings	
	# of reactions
1	37

- ✓ 2. Combine the following reagents to prepare the master mix (contains 0.5 additional reaction to account for pipetting discrepancies):

Master mix		
	A	B
1		37.5 X [µL]
2	ddH2O	1068.75
3	10X LA PCR Buffer	187.5
4	25 mM MgCl2	187.5
5	dNTP Mix (2.5 mM each)	300
6	10 µM OLPr008 MAT -33 UKO F	37.5
7	10 µM OLPr012 MAT +34 SEY UKO R	37.5
8	LA Taq Polymerase (5 U/µL)	18.75
9	OLP045 (pFA6a-His3MX6) (5 ng/µL)	37.5
10	Total	1875

- ✓ 3. PCR program: 94°C/2 min—[94°C/0.5 min—61°C/0.5 min—68°C/1.5 min]x30—68°C/5 min—10°C/∞

Gel


- ✓ 4. Gel: sample volume see below/0.5 µg O'GeneRuler DNA Ladder Mix/1% TAE/120 V/35 min

Load 0.1, 0.2, 0.4 and 0.8 µL of the PCR product (i.e. by diluting PCR product 1:10 by mixing 2 µL product with 18 µL H₂O, then load 1, 2, 4, and 8 µL of the diluted product).

Load 10 µL of no-template control.

Expected product size: 1447 bp.

Estimated concentration of PCR product: 60 ng/0.6 µL = 100 ng/µL.



Additional Resources

Primers and template

	A	B
1	Forward Primer	OLPr008 MAT -33 UKO F
2	Reverse Primer	OLPr012 MAT +34 SEY UKO R
3	Template	OLP045 (pFA6a-His3MX6)

PCR Volume and reagent concentrations

	A	B
1	PCR Volume [μ L]	50
2	Starting buffer conc [X]	10
3	Final buffer conc [X]	1
4	Starting Mg ²⁺ conc [mM]	25
5	Final Mg ²⁺ conc [mM]	2.5
6	Starting dNTP conc each [mM]	2.5
7	Final dNTP conc each [mM]	0.4
8	Starting Primer F conc [μ M]	10
9	Final Primer F conc [μ M]	0.2
10	Starting Primer R conc [μ M]	10
11	Final Primer R conc [μ M]	0.2
12	Starting Polymerase conc [U/ μ L]	5
13	Final Polymerase conc [U/ μ L]	0.05
14	Starting Template conc [ng/ μ L]	5
15	Final Template conc [ng/ μ L]	0.1

Single PCR reaction

	A	B
1		[μ L]
2	ddH ₂ O	28.5
3	10X PCR Buffer	5
4	25 mM MgCl ₂	5
5	dNTP Mix (2.5 mM each)	8
6	10 μ M Primer F	1
7	10 μ M Primer R	1
8	LA Taq Polymerase (5 U/ μ L)	0.5
9	Template (5 ng/ μ L)	1
10	Total	50